


# The activity of fibroblast activation protein (FAP) is reflected by a specific fragment of type III collagen that can be serologically assessed and serve as a non-invasive biomarker

Rasmus Sund Pedersen<sup>1</sup>, Theodora Chrysoulidou<sup>1</sup>, Morten A. Karsdal<sup>1</sup>, Nicholas Willumsen<sup>1</sup>

<sup>1</sup>Nordic Bioscience, Herlev, Denmark,


## Background

- Non-small cell lung cancers (NSCLC) are associated with both inflammatory and fibrotic tumor microenvironments (TME).
- Inflammation is associated with an influx of macrophages, which are releasing matrix metalloproteinases (MMPs).
- MMPs are cleaving extracellular matrix proteins including type III collagen, which results in release of the C3M fragment into circulation.
- Tumor fibrosis is associated with accumulation of fibroblasts and high expression of Fibroblast activation protein (FAP), which is low or absent in benign tissue.
- FAP has unique proteolytic activity and have been shown to cleave various ECM proteins including collagens.
- We hypothesized that FAP-cleaved type III collagen would result in release of a specific fragment (C3F) into circulation, that could reflect FAP-activity in the TME.




## Methods

- To quantify FAP activity, we use *ProteinFingerPrint technology™* to develop a neo-epitope specific CLIA targeting C3F, a fragment specifically generated from FAP-mediated cleavage of type III collagen.
- Specificity of the biomarker towards FAP activity was compared by cleaving MMP13 pre-cleaved type III collagen with FAP or MMP9 followed by measurements of C3F and C3M.
- C3F and C3M were measured in serum from 65 patients with NSCLC and the correlation between the two biomarkers was investigated using Spearman's correlation.
- C3F was also measured in a cohort of 32 patients with NSCLC and 31 healthy individuals and the difference in C3F levels between the two groups was evaluated using Mann-Whitney test and AUROC.




### C3F reflects FAP activity

The C3F fragment could only be detected with FAP-cleaved type III collagen and not with full length or MMP9-cleaved type III collagen.




The C3M fragment could only be detected with MMP9-cleaved type III collagen and not with full length or FAP-cleaved type III collagen.





### No correlation between C3F and C3M was observed

Serum levels of C3M and C3F in patients with NSCLC did not correlate, thus endorsing the claim that the two cleavage-fragments reflect different biological processes.



### C3F showed diagnostic potential for patients with NSCLC

C3F was significantly elevated in serum from patients with NSCLC compared to healthy controls ( $p < 0.0001$ ) and showed diagnostic accuracy of 98%.



## Conclusion

FAP activity can be assessed by targeting a FAP-cleaved fragment of type III collagen, C3F and function as a non-invasive biomarker for patients with NSCLC. By showing that this fragment differs from the MMP-generated fragment C3M, this study supports the notion that different biological processes are reflected in the proteolytic degradation of the ECM and can be reflected by specific circulating protein fragments.

**Contact:** Rasmus Sund Pedersen, [rap@nordicbio.com](mailto:rap@nordicbio.com)