

NordicPRO-C3™

Kit in a box

NordicPRO-C3™ is a highly reliable prognostic biomarker assay that measures active fibrosis by detecting the formation of type III collagen—a key driver of fibrogenesis across major organs. NordicPRO-C3™ can be utilized in the liver, kidney, lungs, heart, intestines, skin, and also in solid tumor research, evidenced by a [Letter of Support by the FDA](#). This biomarker assay can also be implemented in fibroblast cultures modelling such diseases.

NordicPRO-C3™ provides a powerful tool for:

- Monitoring treatment efficacy and pharmacodynamics
- Tracking disease progression over time
- Actionable insights into drug mode of action
- Direct measure of fibroblast activity across fibrotic conditions

With our new ready-to-use [Kit in a box](#), we aim to directly support laboratory scientists and academic researchers in their investigative efforts.



Which sample types are compatible with this kit?

- Human cell culture medium
- Human biological fluids (e.g., serum, EDTA plasma, citrate plasma)
- In vitro cell-culture supernatants (e.g., fibroblast cultures)
- Ex vivo model matrices (e.g., conditioned medium from precision-cut tissue slice cultures)



How many samples can be analyzed with this kit?

- Each 96-well plate accommodates approximately 80 single measurements or 40 duplicate samples. The minimum purchase requirement is two kits, enabling up to 160 single or 80 duplicate measurements.

Order now